

The new alpha trace EMG generation

Software modules

NCS module

Motor and sensory NCS,
Repetitive stimulation, F-wave,
H-reflex, m/s inching,
Blink-reflex, collision

EMG module

Spontaneous, interference Multi-MUP,
Macro-EMG, jitter (single fiber),
MUNE, tremor, galvanic skin response, HRV

EP module

AEP - late, mid-latency and brainstem
SEP - early and late latency
VEP - pattern monitor, flash
MEP- if a magnetic stimulator is available

Contact us for a free demo :
office@alphatrace.at

All-in one PC

- PC and screen all-in one
- Excellent legibility and overview
- Height adjustable

Superlight system cart

- Flexible and open design
- Height adjustable with one finger only
- For sitting or standing operation

Evidence EMG

- Mobile and compact design
- Integrated keyboard with shortcuts
- Specialized and modern software
- Reliable signal quality



EMG accessories catalogue
Get a sample for free:
office@alphatrace.at

Additional equipment

- Tendon hammer
- Foot switch
- Temperature sensor
- Adjustable electro stimulating probe
- Magnetic stimulator
- USB-/Bluetooth dedicated keyboard

Models

Evidence 3102 evo

- 2-channel amplifier
- Only 320 grams light
- Dedicated controls on the device
- Customizable function keys
- Sharp display



Evidence 4004 evo

- 4-channel amplifier
- Expandable to 8 channels
- Ergonomic dedicated keyboard
- Unique confidence analysis for excellent results in evoked potentials



Evidence 5005 evo

- 5-channel amplifier
- Compact design - ideal for mobile use
- Power supply via USB
- Low number of connection cables

